

Fax back Re-order Form

Fax back to Commercial Microbiology Ltd on +44 (0)1224 706012.

Important - Please use the following codes to determine the salinity required for each media:

F/W (Freshwater) : <10,000 TDS **S/W** (Seawater) : 32-34,000 TDS
P/W (Produced Water) : 60,000 TDS **B** (Brine) : >100,000 TDS

<u>ITEM</u>	<u>Salinity Code</u>	<u>Number Required</u>
MPB		
SRB/2		
PRD		
API		
NRB		
NOB		
FAB		
SRB Rapid-Chek Kit	-	
Biocoupons	-	
Biostuds	-	
Plastic Sidestream®	-	
St. Steel Sidestream®	-	
<u>Other (please state below)</u>		
NB – All media comes in packs of 12 vials.		

Signed - _____
Print Name - _____
Company - _____
Telephone Number – _____
Date - _____
Purchase Order Number - _____

Some important information regarding the supply, distribution and use of bacterial growth media produced and supplied by Commercial Microbiology @.

- Media should be stored away from direct sunlight.
- Orders should be checked on arrival, and any shortages / breakages must be reported within 3 days.
- Any correspondence regarding media after its' expiry date will not be entered into.
- Slight differences in colour between batches of the same media is common, and should be expected. Upon extended storage, it is normal for the media colour to change slightly.
- All media undergoes positive release. For special orders that are required urgently, orders may be dispatched before quality checks have been completed. The client will be responsible for the cost of replacement should this later be necessary.
- SRB media containing a nail is produced anaerobically. On occasions, oxygen may enter the vials upon storage, which then reacts with the nail, and is thereby removed from solution. This can cause rusting of the nail, and the media may become discoloured. Such a reaction returns the media to its anaerobic state, and is in no way detrimental to results.